



First, select blocks from your block library by clicking on them. Then, use the Export Settings form to download starter code for selected blocks.

### Block Selector


Select Clear Selected




every



field



bme280



transmit\_msg

### Export Settings

Currently Selected:  
every, field, bme280,  
transmit\_msg

Block Definition(s)  
Format: JSON  
File Name:

Generator Stub(s)  
Language: JavaScript  
File Name:

Export

### Export Preview

Block Definitions:

```
[{
  "type": "every",
  "message0": "every %1 seconds %2 %3",
  "args0": [
    {
      "type": "field_number",
      "name": "DURATION",
      "value": 10,
      "min": 0
    },
    {
      "type": "input_dummy"
    },
    {
      "type": "input_statement",
      "name": "STMTS"
    }
  ],
  "colour": 120,
}
```

Generator Stubs:

```
Blockly.JavaScript['every'] = function(block) {
  var number_duration = block.getFieldValue('DURATION');
  var statements_stmts = Blockly.JavaScript.statementToCode(t
  // TODO: Assemble JavaScript into code variable.
  var code = '...;\n';
  return code;
};

Blockly.JavaScript['field'] = function(block) {
  var text_name = block.getFieldValue('NAME');
  var value_name = Blockly.JavaScript.valueToCode(block, 'nan
  // TODO: Assemble JavaScript into code variable.
  var code = '...';
  // TODO: Change ORDER_NONE to the correct strength.
  return [code, Blockly.JavaScript.ORDER_NONE];
};

Blockly.JavaScript['bme280'] = function(block) {
  var dropdown_field = block.getFieldValue('FIELD');
```

